

1



2



3



2

3

4



5

6

7

8

9

10

Structural and Socket Components - Paediatric Components

1

2

3

4

5

6

7

8

9

10

Paediatric Aluminium 22mm Structural Components

| Part Number | Description | Weight Limit |
|-------------|---|--------------|
| 1 | NJR-7 40L Paediatric 4 Prong Pyramid Adapter | 60kg |
| 2 | NJR-10T2 250mm Length 22mm Pylon | 60kg |
| 3 | NJR-10TC2 250mm Length 22mm Pylon with Titanium Receiver | 60kg |
| | NJR-10TC2-I 250mm Length 22mm Pylon with Short Adapter | 60kg |
| 4 | NJR-10TC2-A Child's 22mm Tube with Adult Pyramid Receiver | 60kg |

Note: Whilst the components showing the blue water symbol are therefore waterproof, please note that these components do not have drainage slots.

1

2

3

4

5

6

7

8

9

10



Structural and Socket Components - Paediatric Components

Weight Limit - 60kg (unless otherwise stated)

Please note that when ordering components from different ranges for the same build, the lowest rated component becomes the maximum weight for the build.

| Paediatric Titanium Pyramid and Receiver Adapters | | | | | |
|---|-------------|--|--------|--------------|------|
| Part Number | Description | Weight | Height | Weight Limit | |
| 1 | NJR-460-TRM | Paediatric Clamping Base with Rotatable Pyramid | 91g | 13mm | 60kg |
| 2 | NJR-465TRF | Paediatric Clamping Base with Rotatable Receiver | 94g | 18mm | 60kg |
| 3 | NJR-460TB | Paediatric 4-Hole Threaded Base with Clamp | 60g | 9mm | 60kg |
| 4 | NJR-00MS | Paediatric Rotatable Pyramid | 31g | 5mm | 60kg |
| 5 | NJR-50FS | Paediatric Receiver, Female (pictured) | 34g | 10mm | 60kg |
| | NJR-50FS-03 | Paediatric Receiver, Female | 39g | 13mm | 60kg |
| | NJR-50FS-06 | Paediatric Receiver, Female | 46g | 16mm | 60kg |
| | NJR-50FS-09 | Paediatric Receiver, Female | 53g | 19mm | 60kg |
| | NJR-50FS-12 | Paediatric Receiver, Female | 59g | 22mm | 60kg |
| 6 | NJR-500P-C | Paediatric SACH Foot Adapter with Single, Centered Anti-Rotation Ridge. <i>Note:</i> Suitable for use with Trulife Child's Play Energy Foot | 58g | 7mm | 60kg |

Note: NJR-460-TRM, NJR-465TRF and NJR-460TB will bolt to NJR-AL-420L (page 259).

Note: Whilst the components showing the blue water symbol are therefore waterproof, please note that these components do not have drainage slots.

1

2

3

4

5

6

7

8

9

10



Structural and Socket Components - Paediatric Adapters

Paediatric Components: Weight Limit - 45kg

Please note that when ordering components from different ranges for the same build, the lowest rated component becomes the maximum weight for the build.

| Socket Adapter with Rotation | | | | |
|------------------------------|-----------|--------|--------------|--|
| Part Number | Material | Weight | Weight Limit | |
| 1 NMA100K4/A | Aluminium | 60g | 45kg | |

| Tube Clamp Adapter with Pyramid | | | | |
|---------------------------------|-----------|--------|--------------|--|
| Part Number | Material | Weight | Weight Limit | |
| 2 NMA110K15/A | Aluminium | 50g | 45kg | |

Paediatric Components: Weight Limit - 60kg

Please note that when ordering components from different ranges for the same build, the lowest rated component becomes the maximum weight for the build.

| Double Adapter - Titanium | | | | | |
|---------------------------|---------------------|--------|--------|--------------|--|
| Part Number | Material | Weight | Height | Weight Limit | |
| 3 NJR-132D | Titanium (pictured) | 62g | 32mm | 60kg | |
| NJR-145D | Titanium | 74g | 45mm | 60kg | |
| NJR-160D | Titanium | 92g | 60mm | 60kg | |

| Child's Tube Clamp with Pyramid Receiver | | | | |
|--|-----------------------|--------|--------------|--|
| Part Number | Material | Weight | Weight Limit | |
| 4 NJR-100C | 22mm Titanium Adapter | 62g | 60kg | |

| Child's Lamination Plate | | | | |
|--------------------------|--|--------|--------------|--|
| Part Number | Material | Weight | Weight Limit | |
| 5 NJR-AL-420L | Aluminium Lamination Plate - supplied with M5 screws. Suitable for use with the Trulife adapter, NAAAK210. | 51g | 60kg | |

| Child's 4-Hole Pyramid Adapter | | | | |
|--------------------------------|---|--------|--------------|--|
| Part Number | Description | Weight | Weight Limit | |
| 6 NJR-400B | Titanium Paediatric 4-Hole Pyramid Adapter | 36g | 60kg | |
| 7 NJR-450BF-M6 | Titanium Paediatric 4-Hole Receiver Adapter | 54g | 60kg | |
| 8 NJR-400A | Titanium Paediatric 4-Hole Pyramid, Axial Offset | 45g | 60kg | |
| 9 NJR-400D | Titanium Paediatric 4-Hole Pyramid, Diagonal Offset | 48g | 60kg | |

| Child's Tube Clamp with Adult Pyramid Receiver | | | | |
|--|--|--------|--------------|--|
| Part Number | Material | Weight | Weight Limit | |
| 10 NJR-100C-A | 22mm Titanium Adapter with 30mm Receiver | 79g | 60kg | |

Note: All paediatric 4 bolt components shown take M5 screws, but several of them are available in M6 screw versions. If required please speak to a Steeper Product Manager for further information.

Note: Whilst the components showing the blue water symbol are therefore waterproof, please note that these components do not have drainage slots.