


Structural and Socket Components - Streifeneder - Cosmetic Foams and Stockings



Cosmetic Soft Foam Cover
For modular transtibial and knee disarticulation prostheses - made from PUR soft foam, applicable for left and right sides and with Ø 19mm bore.

Part Number	Length
1 NMA600K7	50cm

Cosmetic Soft Foam Cover
For modular transfemoral and hip disarticulation prostheses - made from PUR soft foam with 30° knee flexion, a very distinctive pre-shaped calf, matching nylon cosmetic stocking (N200K15) and with Ø 19mm bore.

Part Number	Side	Calf Circumference
2 NMA600K32/L	Left	32cm
NMA600K32/R	Right	32cm

Cosmetic Stockings
For above-knee paediatric prostheses made from 100% polyamide (nylon).

Part Number	Length	Foot Length
3 NMA200K15/1	34cm	10cm
NMA200K15/2	37cm	11cm
NMA200K15/3	44cm	12cm