



Child's Play Lamination Block - Titanium
If lamination is thick for strength, a longer screw may be required.

Part Number	Description	Weight Limit	Weight
1 NSLB240	Lamination Block	65kg	93g

Child's Play Foot Spacers

Part Number	Foot Size	Height	Weight Limit	Weight
2 NSFS240	13-15cm	6mm	31kg	6g
NSFS241	13-15cm	25mm	31kg	20g
NSFS242	16-18cm	6mm	41kg	9g
NSFS243	16-18cm	25mm	41kg	28g
NSFS244	19-21cm	6mm	65kg	11g
NSFS245	19-21cm	30mm	65kg	45g

Child's Play Ankle Block

Part Number	Foot Size	Weight Limit	Weight
3 NSAB314	13-15cm	35kg	123g
NSAB315	16-18cm	45kg	146g
NSAB316	19-21cm	65kg	179g

Child's Play Symes Nut

Part Number	Foot Size	Weight Limit	Weight
4 NSSY301	13-15cm	31kg	14g
NSSY302	16-18cm	40kg	20g
NSSY303	19-21cm	65kg	28g

Child's Play Lamination Core

Part Number	Foot Size	Weight Limit	Weight
5 NSLC301	13-15cm	31kg	28g
NSLC302	16-18cm	40kg	50g
NSLC303	19-21cm	65kg	79g