

Renegade® AT & Renegade® AT LP

The Renegade AT (All Terrain) is an innovative shock-absorbing prosthetic foot providing unparalleled performance in a lightweight, durable design, which is comfortable and versatile, making it suitable for everyday use.

The full-length sole plate and composite Z-Shock heel technology, controls the energy storage and return to enhance the heel dynamics and propel the user through the gait cycle, providing an incredibly smooth gait at any speed and unmatched performance for today's active amputee.

The Renegade AT (All Terrain) Low Profile is the lightest, most responsive and complete shock system for this build height, making it ideal for the active amputee with a longer residual limb. Featuring a full-length sole plate and composite Z-Shock heel technology, the Renegade AT LP prosthetic foot provides an incredibly smooth gait at any speed.

Features and Benefits

- Carbon fibre Z-shaped heel is durable and maintenance free
- Lightweight solutions at 1/3 the weight of traditional shock foot products
- Low profile
- Patented Z-Shock Technology™
- EnduraCore® Technology to improve durability
- Available in 9 stiffness categories



Water resistant - safe for contact with water and submersion*

Cosmetics



Mid profile foot shell



Ordering Information

Trial Period - This foot has a 30 day trial period - only on sizes 23-30 and categories 2-8.

Impact Level Descriptions

Moderate: Daily activities include up to unlimited walking, climbing stairs and occasional moderate lifting (this does not include running).

High: Daily activities include fast walking, jogging, running, lifting heavy objects and/or recreational sports.

Ordering Examples - Prosthetic Foot

Part Number	Category	Size	Left / Right
Renegade AT - NR11526R			
NR11	5	26	R
Renegade AT LP - NR16526R			
NR16	5	26	R

Includes foot shell and Spectra® Sock

Category Selection

Weight (kg)	Impact Level		
	LOW	MOD	HIGH
44-52	1	1	2
53-59	1	2	3
60-68	2	3	4
69-77	3	4	5
78-88	4	5	6
89-100	5	6	7
101-116	6	7	8
117-130	7	8	9
131-147	8	9	
148-166	9		

Technical Information

Build Height and Weight

Size (cm)	Build Height (cm)	Foot inc. Shell
Renegade AT		
22-25	16.2	625g
26-28	18.4	765g
29-31	20	945g
Renegade AT LP		
22-25	12.7	625g
26-28	14	760g
29-31	14.9	930g

Product Specification

Weight Limit	Impact Level	Weight of Unit	Sizes (cm)	Build Height (cm)	Warranty
Renegade AT					
166kg		625g (Size 25)	22-31	16.2 (Size 25)	3 years Foot shell - 6 months
Renegade AT LP					
166kg		625g (Size 25)	22-31	12.7 (Size 25)	3 years Foot shell - 6 months

*Water note: Rinse with fresh water and dry following exposure to water, particularly salt or chlorinated.

Spare Parts Ordering Information

Part Number	Foot Shell Type
NFTC-3M-0▲●4-R■	Standard, with CAP
NFTC-3M-1▲●4-R■	Standard, without CAP
NFTC-2H-1▲●4-R■	High-profile, without CAP

Part Number	Cosmetic attachment plate
NCAP-3M-PO-▲●4-■	

Spare Parts Ordering Guide

- ▲ = Size in cm (22-31)
- = Skin Tone (L = Light, M = Medium, D = Dark)
- = "L" or "R" for Left or Right

Spectra® Sock	Size (cm)	Sock Size (cm)
NSO-NPS-X0030-00	22-25	30
NSO-NPS-X0035-00	26-28	35
NSO-NPS-X0040-00	29-31	40

X denotes colour. Insert 1 for white, or 2 for black. **Note:** Default colour is black

Bumper Kit	Size (cm)
NKIT-00-11600-00	22-25
NKIT-00-11700-00	26-28
NKIT-00-11800-00	29-31