

Seattle Child's Play® Energy Foot

The Child's Play Energy foot is a paediatric version of the Seattle Energy. Its integral cosmesis varies by size to reflect the natural appearance of a child's maturing foot and provides shock absorption at heel strike, whilst the composite keel allows a smooth roll-over with dynamic response at toe-off for a natural gait. Depending on the child's weight and foot size, there are six keel options to choose from to help optimise the foot's action.



***Water-resistant:** Safe for contact with water and occasional submersion. Rinse immediately with fresh water and dry following exposure to water, particularly salt or chlorinated.

Features and Benefits

- Low profile selectable keel design
- Smooth roll-over and dynamic response
- Natural appearance
- Lightweight
- Choice of foot shell colours
- Sandal toe (sizes 16-21)
- Large range of sizes (13-21)
- Safe for occasional contact with water, but not intentional, frequent contact or submersion. After contact with salt or chlorinated water, foot should be rinsed/cleaned.

Cosmetics



Light/dark Sandal Toe Foot shell



Foot construction with foot shell removed



Paediatric - Trulife - Child's Play® Energy Foot

Ordering Information

There is no trial period on this foot.

Part Number	Description	Weight Limit
NSEF132	Child's Play Energy Foot, Light	65kg
NSEF144	Child's Play Energy Foot, Dark	65kg

Ordering Example - NSEF132-18-R-4

Part Number	Size	Left / Right	Category
NSEF132	18	R	4

Spare Parts Ordering Information

Part Number	Description	Weight Limit
NSFB546	Bolt Kit Childs Play Energy, 8mm x 80mm	65kg
NSSY301	Child's Play Symes Nut, 13-15cm	31kg
NSSY302	Child's Play Symes Nut, 16-18cm	40kg
NSSY303	Child's Play Symes Nut, 19-21cm	65kg

Category Selection Instructions

To determine which keel is appropriate for your patient, check the below table for keel availability by foot size. Please note that size 13-15 feet do not have sandal toe foot shells.

Important Note: Choose the next higher keel level if your patient has a long BK limb, regularly carries heavy loads, or is at the upper limit for weight or impact level.

Warning: Choosing a lower strength category than what is suggested based on the above procedure and patient data will void the warranty. If your patient's weight exceeds the limits of the chart, please browse the catalogue for other foot options or contact our customer services department.

Weight Limit	Keel	Foot Size									
		13	14	15	16	17	18	19	20	21	
<15kg	1	✓	✓	✓	-	-	-	-	-	-	-
16-25kg	2	✓	✓	✓	✓	✓	✓	-	-	-	-
26-35kg	3	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
36-45kg	4	-	-	-	✓	✓	✓	✓	✓	✓	✓
46-55kg	5	-	-	-	-	-	-	✓	✓	✓	✓
56-65kg	6	-	-	-	-	-	-	✓	✓	✓	✓

Product Specification

Weight Limit	Impact Level	Weight	Sizes (cm)	Build Height (cm)	Warranty
65kg	L M H	150g Size 17	13-21	4.5 Size 17	2 years Foot shell - 1 year
Size	Minimum Clearance (cm)		Heel Rise (cm)		
13-15	4		0.63		
16-18	4.5		0.95		
19-21	4.5		1.27		