

Seattle Natural

The Seattle Natural foot, with its Delrin Keel, has a slightly more responsive action than the SACH, making it suitable for low to medium-low impact levels, single speed, home ambulators or new amputees (K1-K2). It is supplied already assembled with M10 pyramid and foot bolt, with a Lightfoot cosmesis, in two colour options.

Features and Benefits

- Lightweight
- One standard Delrin keel per foot size
- Propels the user through cadence
- Natural looking Lightfoot cosmesis
- Available in both Light and Dark foot cosmesis
- Sandal toe only
- Safe for occasional contact with water, but not intentional, frequent contact or submersion. After contact with salt or chlorinated water, feet should be rinsed/cleaned.

Cosmetics



Ordering Information

Trial Period - There is no trial period on this foot.

Impact Level Descriptions

Low: Low speed walking and low intensity activities - no variable cadence.

Part Number	Description	Available Sizes (cm)
NSNF160	9.5mm Heel Light, Lightfoot cosmesis	22-30 cm
NSNF163	9.5mm Heel Dark, Lightfoot cosmesis	22-30 cm

Ordering Example - NSNF160-25L


Part Number	Size	Left / Right
NSNF160	25	L

Feet - Trulife - Seattle Natural

Product Specification

Build Height and Heel Rise			
Part Number	Size (cm)	Build Height (cm)	Heel Rise (mm)
NSNF160	22-30	8.5	9.5
NSNF163	22-30	8.5	9.5

Weight Limit	Impact Level	Weight of Unit	Sizes (cm)	Build Height (cm)	Warranty
102kg		448g Size 25	22-30	8.5	2 years Foot shell - 1 year



Seattle Natural keel view