

Seattle SACH

The Seattle SACH foot is suitable for new amputees who are likely to ambulate with low to medium-low impact levels (K1). It is supplied already assembled with M10 pyramid and foot bolt, with a Lightfoot cosmesis in two colour options.

Features and Benefits

- Easy roll-over
- Lightweight
- Natural looking Lightfoot cosmesis
- Sandal toe
- Available in both Light and Dark foot cosmeses
- Safe for occasional contact with water, but not intentional, frequent contact or submersion. After contact with salt or chlorinated water, feet should be rinsed/cleaned.

Cosmetics



Ordering Information

Trial Period - There is no trial period on this foot.

Impact Level Descriptions

Low: Low speed walking and low intensity activities - no variable cadence.

Part Number	Description	Available Sizes (cm)
NSCH110	9.5mm Heel Light, Lightfoot cosmesis	22-30 cm
NSCH113	9.5mm Heel Dark, Lightfoot cosmesis	22-30 cm

Ordering Example - NSCH110-25L


Part Number	Size	Left / Right
NSCH110	25	L

Feet - Trulife - Seattle SACH

Product Specification

Build Height and Heel Rise			
Part Number	Size (cm)	Build Height (cm)	Heel Rise (mm)
NSCH110	22-30	8.5	9.5
NSCH113	22-30	8.5	9.5

Weight Limit	Impact Level	Weight of Unit	Sizes (cm)	Build Height (cm)	Warranty
102kg		549g Size 25	22-30	9 Size 22-26	2 years Foot shell - 1 year



Seattle SACH keel view