

Senator®

The Senator is a lightweight, economical, energy returning foot that shares composite engineering features with the advanced Freedom series. It is designed to meet the needs of low to moderately active amputees seeking a prosthetic foot that will allow them to ambulate at varied speeds and participate in daily activities such as cycling, golfing, hiking, or even light jogging. The Senator is a durable foot that is simple to fit and finish.

CES FRE FO1: Use this code to read the Senator clinical evaluation summary on the Steeper website.

Features and Benefits

- Lightweight, energy-return design
- Integral male pyramid
- Low maintenance, durable design
- Available in 6 stiffness categories



Water resistant - safe for contact with water and submersion*

Cosmetics



Mid profile foot shell



Ordering Information

Trial Period - There is no trial period on this foot.

Impact Level Descriptions

Low: Daily activities include mostly level ground walking, moving around in the home and the community.

Moderate: Daily activities include up to unlimited walking, climbing stairs and occasional moderate lifting (this does not include running).

Ordering Example - NVS1324L

Part Number	Category	Size	Left / Right
NVS1	3	24	L

Includes foot shell and Spectra® Sock


Category Selection

Weight (kg)	Impact Level	
	LOW	MOD
45-55	1	2
56-66	2	3
67-80	3	4
81-96	4	5
97-114	5	6
115-136	6	

Build Height and Weight

Size (cm)	Build Height (cm)	Foot inc. Shell
22-25	13	520g
26-28	14	655g
29-31	15.2	780g

Technical Information



Product Specification

Weight Limit	Impact Level	Weight of Unit	Sizes (cm)	Build Height (cm)	Warranty
136kg	L M H	520g	22-31	13 (size 25)	2 years Foot shell - 6 months

*Water note: Rinse with fresh water and dry following exposure to water, particularly salt or chlorinated.

Spare Parts Ordering Information



Part Number	Foot Shell Type
NFTC-3M-0▲●4-R■	Standard, with CAP
NFTC-3M-1▲●4-R■	Standard, without CAP
NFTC-2H-1▲●4-R■	High-profile, without CAP



Cosmetic Attachment Plate
NCAP-3M-PO-▲●4-R■

Spare Parts Ordering Guide

- ▲ = Size in cm (22-31)
- = Skin Tone (L = Light, M = Medium, D = Dark)
- = "L" or "R" for Left or Right

Spectra® Sock	Size (cm)	Sock Size (cm)
NSO-NPS-X0025-00	22-23	25
NSO-NPS-X0030-00	24-26	30
NSO-NPS-X0035-00	27-28	35
NSO-NPS-X0040-00	29-31	40

X denotes colour. Insert 1 for white, or 2 for black. **Note:** Default colour is black

Bumper Kit	Size (cm)
NKIT-00-11383-00	22-25
NKIT-00-11384-00	26-28
NKIT-00-11385-00	29-31