

Velocity®

Delivering a lower profile than the Soleus, the recently redesigned and reinforced Velocity features coupled toe springs that work together to provide a similarly progressive and smooth roll-over.

Velocity combines high functionality and low maintenance, with ease of alignment and set up. The additional inserts provided allow fine tuning of the foot function.

Features and Benefits

- Integrated Spring Technology (IST)
- Multi-axial function
- IntelliWeave composite to bring maximum strength and flexibility
- Easily adjustable tuning with three different heel wedges and a toe wedge provided

Cosmetics



CAP



Foot shell is anatomically designed



Ordering Information

Trial Period - This foot has a 30 day trial period.

Impact Level Descriptions


Low: Daily activities include mostly level ground walking, moving around in the home and the community.

Moderate: Daily activities include up to unlimited walking, climbing stairs and occasional moderate lifting (this does not include running).

High: Daily activities include fast walking, jogging, running, lifting heavy objects and/or recreational sports.

Velocity foot includes sock and foot shell. See Order Forms on page 38 for additional ordering information.

Technical Information



Product Specification

Mounting	Weight Limit	Impact Level	Weight of Unit	Sizes (cm)	Build Height (cm)	Warranty	Foot Shell
Endo (IP)	113kg (Sizes 25-26) 125kg (Sizes 27-30)		631g (Size 26)	25-30	12 (Sizes 25-26) 13.2 (Sizes 27-30)	3 years Foot shell - 6 months	 Caucasian, Tan, Brown or Jet Black Narrow or wide width

Shelltread - Ordering Information

Shelltread	Size
NCST	25-30cm

Shelltread Ordering Guide	
	= "L" or "R" for Left or Right
	= Colour (N = Neutral, J = Jet Black)
	= Foot size in cm (25-30)
	= Wide

Feet - College Park - Velocity

Spare Parts Ordering Information



Foot Shell	Size
NV S	25-30
NV S	25-30

W = Wide



CPI Sock	Size
NV SS 2526	25-26
NV SS 2730	27-30

Spare Parts Ordering Guide

- = "L" or "R" for Left or Right
- = Colour (C = Caucasian, T = Tan, B = Brown, J = Jet Black)
- = Foot size in cm (25-30)
- = Wide