

The Zumo foot is designed to withstand 205kg, making it ideal for heavier users.

The split composite keel provides increased dynamic response, stability and improved medial/lateral compliance. The removable sandal toe foot shell, which comes in two colours, offers a high quality, natural cosmesis.

Features and Benefits

- Heavy-duty foot withstands up to 205kg of weight
- Durable composite design for low to moderate users
- Pre-assembled with a male pyramid for quick and easy installation
- Split composite keel for increased stability
- Dynamic toe and cushioned heel response
- Removable, lifelike, sandal toe foot shell available in light or dark
- Safe for occasional contact with water, but not intentional, frequent contact or submersion. After contact with salt or chlorinated water, the foot should be rinsed/cleaned

Cosmetics



Light/dark Sandal Toe Footshell



Foot construction with foot shell removed



Ordering Information

Trial Period - There is no trial period on this foot.

Impact Level Descriptions

Low: Low speed walking and low intensity activities - no variable cadence.

Medium: Medium impact activities including light sports.

Part Number	Description	Available Sizes (cm)
NSZM450	9.5mm Heel, Light	25-30
NSZM453	9.5mm Heel, Dark	25-30

Ordering Example - NSZM450-25-L

Part Number	Size	Left/Right
NSZM450	25	L
NSZM453	25	L

Spare Parts Ordering Information

Foot Shell Ordering Example - NSFC19025L

Part Number	Size	Left/Right
NSFC190	25	L

Replacement Foot Shell, Light		
Part Number	Size (cm)	
NSFC190	25-30	

Replacement Foot Shell, Dark		
Part Number	Size (cm)	
NSFC193	25-30	

Product Specification

Weight Limit	Impact Level	Weight	Sizes (cm)	Build Height (cm)	Warranty
205kg		750g Size 25cm	25-30	9 Size 25cm	2 years Foot shell - 6 months

Compatible Components



Heavy Duty Tube Clamp with 30mm Pyramid Receiver		
Part Number	Weight Limit	Man.
NHD-100C-A	227kg	APC



Tube Clamp Adapter (34mm)		
Part Number	Weight Limit	Man.
NMA11A225/T	227kg	SF



Light Aluminium Alloy Tube (34mm)		
Part Number	Weight Limit	Man.
NMA28A3/A300	227kg	SF



Laminating Socket Adapter		
Part Number	Weight Limit	Man.
NMA10A7/T	227kg	SF

Manufacturer:

APC = American Prosthetic Components

SF= Streifeneder

Additional Parts Ordering Information

For more heavy duty components please refer to page 260.