

The College Park Genesee liner offers all the benefits of a gel liner in an affordable package. This durable liner is shaped for fit and flexibility, with cushioning comfort and excellent limb conforming properties. It fits a wide range of users.

The unique fabric improves the fibre-to-gel bond, reduces pistoning, and makes donning the socket easier. Mineral oil is infused into the liner to help absorb shear forces and maintain moisture.

Features and Benefits

- Retro-fits with industry standard liners of equal thickness and profile
- Infused with mineral oil to moisturise the skin and absorb shear forces
- Unique fabric provides a superior fibre to gel bond, reducing pistoning and easing socket donning
- 41cm long for the perfect fit
- Activity Level K1-K3

Ordering Information


Ordering Example: NGLK-G-6C-TM

Part ID	Thickness	Type	Profile	Size
NGLK-G				
	3mm 3	Pinlock L	Tapered T	Medium M
	6mm 6	Cushion C	Uniform U	Medium Plus MP
	9mm 9			Large L
				Large Plus LP
				Extra Large XL

Size Information

Thickness	3mm		6mm		9mm		Proximal Circ.	Distal Circ.
	Uniform	Tapered	Uniform	Tapered	Uniform	Tapered		
Medium	✓	✓	✓	✓	✓	X	23 - 33cm	20 - 28cm
Medium+	✓	✓	✓	✓	✓	X	33 - 43cm	20 - 28cm
Large	✓	X	✓	✓	✓	X	36 - 46cm	23 - 30cm
Large+	✓	X	✓	X	✓	X	41 - 53cm	28 - 38cm
Extra Large	✓	X	✓	X	X	X	43 - 69cm	30 - 51cm

✓ Available X Not available



Technical Information

Material	Thermoplastic Elastomer (TPE) Gel
Length	40.6cm (+/- 1.2cm)
Type	Cushion or Locking
Size	M, M+, L, L+, XL
Profile	Tapered or Uniform
Thickness	3mm, 6mm, 9mm
Activity Level	K1-K3
Warranty	6 months

