

Genesee Liner

The College Park Genesee liner provides the ideal benefits required of a gel liner in an affordable package. Shaped for fit and flexibility with cushioning comfort and excellent limb conforming properties, this durable liner accommodates a wide range of users.

The unique fabric provides a superior fibre to gel bond, reduces pistoning and eases donning. The liner is infused with mineral oil to help absorb shear forces and provide continuous moisture.

Features and Benefits

- Retro-fits with industry standard liners of like thickness and profile
- Infused with mineral oil to moisturise the skin and absorb shear forces
- Unique fabric provides a superior fibre to gel bond, reducing pistoning and easing donning
- 41cm long for the perfect fit
- Activity level K1-K3

Ordering Information

Part Number

NGLK-G ■●-◆▲

- - Thickness (3 = 3mm, 6 = 6mm, 9 = 9mm)
- - Type (L = Pinlock, C = Cushion)
- ◆ - Profile (T = Tapered, U = Uniform)
- ▲ - Size (M = Medium, MP = Medium Plus, L = Large, LP = Large Plus, XL = Extra Large)

Technical Information

Material	Thermoplastic Elastomer (TPE) Gel
Length	40.6cm (+/- 1.2cm)
Type	Cushion or Locking
Size	M, M+, L, L+, XL
Profile	Tapered or Uniform
Thickness	3mm, 6mm, 9mm
Warranty	6 months

Size Information

Thickness	3mm		6mm		9mm		Proximal Circ.	Distal Circ.
	Uniform	Tapered	Uniform	Tapered	Uniform	Tapered		
Medium	✓	✓	✓	✓	✓	X	23 - 33cm	20 - 28cm
Medium+	✓	✓	✓	✓	✓	X	33 - 43cm	20 - 28cm
Large	✓	X	✓	✓	✓	X	36 - 46cm	23 - 30cm
Large+	✓	X	✓	X	✓	X	41 - 53cm	28 - 38cm
Extra Large	✓	X	✓	X	X	X	43 - 69cm	30 - 51cm

✓ Available X Not available

