

Uniprox Liners

SoftSkin and SoftSkin Air

SoftSkin S30	SoftSkin Air S30	SoftSkin S40	SoftSkin Air S40	SoftSkin S50	SoftSkin Air S50	SoftSkin TF	SoftSkin Air TF
The gentle liner	The breathable silicone-liner	The comfortable liner	The breathable silicone-liner	The active liner	The breathable silicone-liner	The comfortable liner	The breathable silicone-liner
							

Shore hardness	30	30	40	40	50	50	40	40
Thickness	3mm	3mm	3mm	3mm	3mm	3mm	3mm	3mm
Length	370mm	370mm	370mm	370mm	370mm	370mm	370mm	370mm
Activity level	1,2,3	1,2,3	1,2,3	1,2,3	1,2,3	1,2,3	1,2,3 and 4	1,2,3 and 4
Reduced perspiration due to micropores		●		●		●		●
Graduated textile cover provides secure distal containment	●	●	●	●	●	●		
Stainless steel pin fitting	●	●	●	●	●	●		
Soft cushioning of the residual limb	●	●	●	●			●	●
Good adherence from the silicone	●		●					
Order Number (With Distal Connector)	NS30-■	NS30-■-AIR	NS40-■	NS40-■-AIR	NS50-■	NS50-■-AIR	NSTF-■	NSTF-■-AIR
Cushion Liner Order Number (Without Distal Connector)	NS30CL-■	NS30CL-■-AIR	NS40CL-■	NS40CL-■-AIR	NS350CL-■	NS50CL-■-AIR		