

SureFit Sleeve



30° Pre-Flexion



Innovative Weave Design



High-Density Gel





The **SureFit Sleeve** is constructed using High-Density TPE Gel (HD Gel) to provide maximum comfort for the patient and long-lasting security. With a 30° pre-flexion it allows for ease when bending and reduces bunching behind the knee; whilst the intelligent weave is adapted to differing areas of stretch and wear.

The SureFit Sleeve HD gel is characterised by a firm texture, though it is more compliant under compression in comparison to silicone. This extreme versatility allows active amputees a greater degree of control and feedback during everyday activities, without creating shear forces across the skin.

Features:

- 30° pre-flexion
- 3mm of uniform thickness
- Low modulus of elasticity for correct compression and stability
- One length (400mm) to accommodate most users in three sizes
- Two colour options available: Grey or Beige

Ordering Information

Circumference	Sleeve Size	Colour	Product Part No.
20 - 34cm	Small	Beige	NSFX20-SM-B
20 - 34cm	Small	Grey	NSFX20-SM-G
31 - 40cm	Medium	Beige	NSFX26-MD-B
31 - 40cm	Medium	Grey	NSFX26-MD-G
35 - 60cm	Large	Beige	NSFX32-LG-B
35 - 60cm	Large	Grey	NSFX32-LG-G

Example Order*

Circumference measurement = 31cm, in Beige. Select size **Small**, product number **NSFX20-SM-B**.

Circumference measurement = 40cm, in Grey. Select size **Medium**, product number **NSFX26-MD-G**.

*Measure the thigh circumference 10cm proximal to the patella. If the measured circumference falls within two sleeve sizes, please order the smaller size.

To place an order, contact customerservices@steepergroup.com or call **+44 (0) 870 240 4133**. For more information about SureFit Sleeves, please contact your local Product Manager.

Unit 3, Stourton Link
 Intermezzo Drive, Leeds, LS10 1DF
 Tel: +44 (0) 870 240 4133
 Email: customerservices@steepergroup.com