

Pro Armour PLUS

The latest product from Pro Armour, the Pro Armour PLUS provides the same first-class protection as the Lite, giving amputees the ultimate freedom to walk footwear-free.

Thanks to the low profile rubberised sole, users can enjoy more footwear-free activities such as walking along the shoreline of the beach or being poolside. The sock not only provides grip, but the neoprene properties ensure components are protected from ingress of dirt.

Features and Benefits

- Splash, sand, dust and dirt proof
- 4mm sole thickness
- Promotes physical activity and positive wellbeing
- Low maintenance - washable & reusable
- Available in multiple sizes
- Suitable for users of all activity levels

Foot Compatibility

The Pro Armour PLUS is generally suitable for most SACH, uniaxial, multiaxial, blades, energy storing and patient adjustable feet.



The Pro Armour PLUS can give users more flexibility and independence to do what they love most



The rubberised, non-slip sole gives users confidence in knowing they're safe from slipping on wet and uneven surfaces

Socks - Pro Armour® - Pro Armour PLUS

Ordering Information

Impact Level Descriptions

Low: Low impact activities include walking on uneven surfaces.

Medium: Medium impact activities including light sports.

High: High impact activities such as running and basketball.

*Available in both left (L) and right (R) variants - please choose **L** or **R** when ordering.

**Foot shell sizing varies - use only this as a guide when placing an order.

Pro Armour PLUS		
Part Number*	Foot Shell Size**	Weight
NPAP22	22	155g
NPAP23	23	155g
NPAP24	24	165g
NPAP25	25	180g
NPAP26	26	180g
NPAP27	27	195g
NPAP28	28	200g
NPAP29	29	220g

Technical Information

	Foot Shell Size	Weight	Colour	Material	Activity Level
NPAP22	22	155g	Black	Neoprene	LMH
NPAP23	23	155g	Black	Neoprene	LMH
NPAP24	24	165g	Black	Neoprene	LMH
NPAP25	25	180g	Black	Neoprene	LMH
NPAP26	26	180g	Black	Neoprene	LMH
NPAP27	27	195g	Black	Neoprene	LMH
NPAP28	28	200g	Black	Neoprene	LMH
NPAP29	29	200g	Black	Neoprene	LMH

