

The Relax range from Össur provides a pain-relieving effect from phantom limb pain, phantom limb sensation and idiopathic residual limb pain. This is due to the metallurgical knitted fabric screening against electromagnetic influences. The Relax Night Care socks are the ideal complement to the Relax liners because they can be worn at night as well.

Umbrellan fabric is a unique, patented knitted fabric with an excellent electromagnetic screen. The materials used on the metallurgical base are hypoallergenic, climate-regulating and have antibacterial properties.


The socks are elasticated with a soft top band providing a snug and comfortable non-slip fit. The Clima-Comfort® technology ensures a high degree of comfort by transporting moisture rapidly from the skin to the surface of the stocking, reducing sweating. At the same time, the skin has an optimal oxygen supply to prevent it drying out.



TEST WITHOUT RISK  
"MONEY BACK  
GUARANTEE"


## Features and Benefits

- Pain-relieving effect due to integral electromagnetic shield
- Clima-Comfort technology - warming and cooling function depending on environmental temperature
- Elastic knit with soft top bands for comfortable, reliable fit without slipping or constricting
- Suitable for wearing day and night or whenever the prosthesis is not being worn



## Ordering Information

Circumference A (Level of top band end of residual limb)	Length	Size	Part Number
20-26cm	35cm	S	N188000S
26-33cm	35cm	M	N188000M
33-42cm	35cm	L	N188000L
42-52cm	35cm	XL	N18800XL
52-62cm	35cm	XXL	N1880XXL



### Size Selection

1. Depending on the desired length (L), measure the circumference (A) 35cm above the end of the residual limb.

2. Residual limb length L and the circumference are the basis for selecting the correct Össur RELAX Night Care size.

## Important information

The Össur Relax Night Care Sock is available in one standard length (35cm). The sock can be carefully sewn over for residual limbs < 30cm.

## Money Back Guarantee

### For Relax Night Care Sock

To help assess if the Relax Umbrellan system will reduce or eradicate your patient's phantom limb pain, we are offering a risk free "Money Back Guarantee" for the Össur Relax Night Care Sock. If after 5 days your patient feels no benefit from the product, you can return it to us for a full refund.

If your patient feels a reduction in pain from the sock, they will receive the same benefit from the liner. Likewise, if no benefit is felt from the sock, a reduction in pain would not be expected from the liner.

## Technical Information

Material	Umbrellan with Clima-Comfort
Size / Length	See ordering information
Amputation Level	Transtibial and knee disarticulation
Cleaning	30° hand wash
Items Included	Umbrellan user set (Pain Diary, Pain Scale)
Colour	Actual colour may vary from the image and from one batch to another