

## Skin-Fix Donning Sleeve

Skin-Fix is to be donned on the residual limb before slipping into the socket. Due to the low friction design, donning of the socket becomes an easy and comfortable operation.

Then, the Skin-Fix can be pulled out through the valve easily.

### Features and Benefits

- Low friction
- Long-life durability
- Simple application

## Socks - Uniprox - Donning Sleeve

### Ordering Information

| Part Number | Size     | Upper Circumference | Lower Circumference | Length |
|-------------|----------|---------------------|---------------------|--------|
| NF011748    | XX-Large | 820mm               | 320mm               | 1200mm |
| NF011749    | X-Large  | 820mm               | 320mm               | 1000mm |
| NF011750    | Large    | 720mm               | 310mm               | 1000mm |
| NF011751    | Medium   | 600mm               | 260mm               | 1000mm |
| NF011752    | Small    | 460mm               | 260mm               | 1000mm |
| NF011753    | X-Small  | 420mm               | 240mm               | 1000mm |

