

Transtibial Compression Socks

Intelligent and effective products with deep functional effects for optimum comfort. The Össur compression seal identifies selected high-tech products from Össur in medicine, lifestyle and sports.

The manufacturing process provides accurately dosed compression on the body - for positive effects on health, performance and well-being.


Features and Benefits

- Soft, smooth material for comfortable all-day wearing
- Silicone beaded BK top-band for superior suspension
- Smaller BK size



Size Selection

Below-knee



Ordering Information

Comfort SCS BK

| Length | Size | Cir. cC (cm) | Cir. cE (cm) | Cir. cF (cm) | Part No. |
|-----------------|------|--------------|--------------|--------------|------------|
| Short 28-34 cm | 0 | 26-28 | 28-30 | 38-40 | N1881 ▲ 10 |
| | 1 | 28-30 | 30-32 | 40-42 | N1881 ▲ 11 |
| | 2 | 30-33 | 32-35 | 42-45 | N1881 ▲ 12 |
| | 3 | 33-36 | 35-38 | 45-48 | N1881 ▲ 13 |
| | 4 | 36-39 | 38-41 | 48-51 | N1881 ▲ 14 |
| | 5 | 39-42 | 41-44 | 51-54 | N1881 ▲ 15 |
| Medium 35-41 cm | 6 | 42-45 | 44-47 | 54-57 | N1881 ▲ 16 |
| | 0 | 26-28 | 28-30 | 38-40 | N1881 ▲ 20 |
| | 1 | 28-30 | 30-32 | 40-42 | N1881 ▲ 21 |
| | 2 | 30-33 | 32-35 | 42-45 | N1881 ▲ 22 |
| | 3 | 33-36 | 35-38 | 45-48 | N1881 ▲ 23 |
| | 4 | 36-39 | 38-41 | 48-51 | N1881 ▲ 24 |
| Long 42-48 cm | 5 | 39-42 | 41-44 | 51-54 | N1881 ▲ 25 |
| | 6 | 42-45 | 44-47 | 54-57 | N1881 ▲ 26 |
| | 0 | 26-28 | 28-30 | 38-40 | N1881 ▲ 30 |
| | 1 | 28-30 | 30-32 | 40-42 | N1881 ▲ 31 |
| | 2 | 30-33 | 32-35 | 42-45 | N1881 ▲ 32 |
| | 3 | 33-36 | 35-38 | 45-48 | N1881 ▲ 33 |
| | 4 | 36-39 | 38-41 | 48-51 | N1881 ▲ 34 |
| | 5 | 39-42 | 41-44 | 51-54 | N1881 ▲ 35 |
| | 6 | 42-45 | 44-47 | 54-57 | N1881 ▲ 36 |

▲ - Denotes compression, add 1 for CCL1 and 2 for CCL2

Technical Information

Indications

- Reduction of oedema/swelling of the lower limb after amputation
- Maintenance of limb volume when the prosthesis is removed
- For fixing soft tissues

Contraindications

- If the medical findings forbid compression (e.g. heart failure)

| | |
|-------------------|--------------------------------|
| Material | Umbrellan® with Clima-Comfort® |
| Size / Length | See Ordering Information |
| Amputation Level | Transtibial |
| Cleaning | At 30° hand wash |
| Compression Class | CCL1 and CCL2 |

Transfemoral Compression Socks

Intelligent and effective products with deep functional effects. The Össur compression seal identifies selected high-tech products from Össur in medicine, lifestyle and sports.

The manufacturing process provides accurately dosed compression on the body - for positive effects on health, performance and well-being.

Features and Benefits

- Soft, smooth material with a comfortable feel for all-day wearing
- Movable AK waistbelt



Clima Comfort




Clima Fresh



Size Selection

Above-knee



Ordering Information

Comfort SCS AK

| Length | Size | Cir. cF (cm) | Cir. cG (cm) | Part No. |
|-----------------|------|--------------|--------------|------------|
| Short 18-22 cm | 1 | 30-32 | 42-45 | N1882 ▲ 11 |
| | 2 | 32-35 | 45-49 | N1882 ▲ 12 |
| | 3 | 35-38 | 49-53 | N1882 ▲ 13 |
| | 4 | 38-41 | 53-57 | N1882 ▲ 14 |
| | 5 | 41-44 | 57-61 | N1882 ▲ 15 |
| | 6 | 44-47 | 61-65 | N1882 ▲ 16 |
| Medium 23-27 cm | 1 | 30-32 | 42-45 | N1882 ▲ 21 |
| | 2 | 32-35 | 45-49 | N1882 ▲ 22 |
| | 3 | 35-38 | 49-53 | N1882 ▲ 23 |
| | 4 | 38-41 | 53-57 | N1882 ▲ 24 |
| | 5 | 41-44 | 57-61 | N1882 ▲ 25 |
| | 6 | 44-47 | 61-65 | N1882 ▲ 26 |
| Long 28-32 cm | 1 | 30-32 | 42-45 | N1882 ▲ 31 |
| | 2 | 32-35 | 45-49 | N1882 ▲ 32 |
| | 3 | 35-38 | 49-53 | N1882 ▲ 33 |
| | 4 | 38-41 | 53-57 | N1882 ▲ 34 |
| | 5 | 41-44 | 57-61 | N1882 ▲ 35 |
| | 6 | 44-47 | 61-65 | N1882 ▲ 36 |
| X-Long 33-37 cm | 1 | 30-32 | 42-45 | N1882 ▲ 41 |
| | 2 | 32-35 | 45-49 | N1882 ▲ 42 |
| | 3 | 35-38 | 49-53 | N1882 ▲ 43 |
| | 4 | 38-41 | 53-57 | N1882 ▲ 44 |
| | 5 | 41-44 | 57-61 | N1882 ▲ 45 |
| | 6 | 44-47 | 61-65 | N1882 ▲ 46 |

▲ - Denotes Compression, add 1 for CCL1 and 2 for CCL2

Socks - Össur - Transfemoral Compression Sock

Technical Information

Indications

- Reduction of oedema/swelling of the lower limb after amputation
- Maintenance of limb volume when the prosthesis is removed
- For fixing soft tissues

Contraindications

- If the medical findings forbid compression (e.g. heart failure)

| | |
|-------------------|-------------------------------------|
| Material | Umbrellan® with Clima-Comfort® |
| Size / Length | See Ordering Information |
| Amputation Level | Knee disarticulation / transfemoral |
| Cleaning | At 30° hand wash |
| Compression Class | CCL1 and CCL2 |