

Hood Socks

Hood Classic Wool Socks

- 100% Merino wool
- Available in 2, 3, 5 or 6 ply
- Made-to-measure available
- Lengths: 8"-34" (20-86cm)
- Diameters: 6"-12" (15-30cm)



Wool

Hood Nylon Sheaths

- 100% Nylon
- Single ply thickness
- 11 sizes available
- Pocket seam
- Wicks away perspiration



Nylon

Ordering Information

Hood Nylon Sheaths - Pocket Seam

Part Number	Description	Size / Length	Top Circ (cm)	Length (cm)
N54011U02	Nylon Sheath	Extra Narrow - Short	22-26	30-35
N54012U02	Nylon Sheath	Narrow - Short	30-36	30-35
N54012U03	Nylon Sheath	Narrow - Regular	30-36	40-46
N54012U04	Nylon Sheath	Narrow - Long	30-36	51-61
N54013U02	Nylon Sheath	Regular - Short	40-46	30-35
N54013U03	Nylon Sheath	Regular - Regular	40-46	40-46
N54013U04	Nylon Sheath	Regular - Long	40-46	51-61
N54013U05	Nylon Sheath	Regular - Extra Long	36-46	66-76
N54014U02	Nylon Sheath	Wide - Short	50-60	30-35
N54014U03	Nylon Sheath	Wide - Regular	50-60	40-46
N54014U04	Nylon Sheath	Wide - Long	50-60	51-61

Hood Classic Wool Socks


Ordering Example - Wool Hood Socks

Part Number	Description
N510 ● ▲ U ●	Wool Sock

● = Ply (2,3,5,6)
 ▲ = Proximal Sock Width in Inches (6"-12")
 ● = Sock length in Inches (8"-34")*
 *Length in increments of 2"

Made to Measure Hood Socks (Wool Only)

Ordering Information	
Patient Name	Address
Prosthetist	
Branch Name	Post Code
Account Number	Comments
Width at Top (inches)	
Length (inches)	
Number of Ply	
Quantity	



Below Knee / Transtibial

Above Knee / Transfemoral

— A Proximal Flat Width of Sock in Inches (Circumference divided by 2)

— B Toe Width of Sock in Inches (Circumference 2" proximal of the distal end divided by 2)

— C Sock Length in Inches (Including roll-over)

Note: Please mark any deviations on the drawing.