

Relax Compression Socks


These compression socks within the Relax range are particularly beneficial for patients with sensitive residual limb conditions. With Umbrellan fabric, these socks are designed to shield against electro-magnetic influences as one of the possible causes of phantom pain.

The Clima Comfort within the knit used has a warming and cooling temperature, depending on the ambient temperature, whilst the highly elastic knitted fabric has soft cuffs for a comfortable fit.

Electromagnetic shielding: As evidenced in Fisher K, Oliver S, Sedki I, Hanspal R. The effect of electromagnetic shielding on phantom limb pain: A placebo-controlled double-blind crossover trial. **Prosthetics and Orthotics International**. 2015.


Features and Benefits

- Soft, smooth material for comfortable all-day wearing
- Comfort zone at the distal end protects the scarred area against cold, knocks, friction
- Secure hold due to non-slip topband




Ordering Information

Below-knee



Length	Size	Cir. cC (cm)	Cir. cE (cm)	Cir. cF (cm)	Part No.
BK 30-38cm	2	28-31	30-35	36-43	N1880422
	3	31-34	33-38	40-47	N1880423
	4	34-37	36-42	44-51	N1880424
	5	37-40	40-45	49-56	N1880425
BK 38-46cm	6	40-44	43-48	53-60	N1880426
	2	28-31	30-35	36-43	N1880432
	3	31-34	33-38	40-47	N1880433
	4	34-37	36-42	44-51	N1880434
BK 38-46cm	5	37-40	40-45	49-56	N1880435
	6	40-44	43-48	53-60	N1880436

Above-knee



Length	Size	Cir. cF (cm)	Cir. cG (cm)	Part No.
AK 20-26cm	2	30-34	42-50	N1880522
	3	34-38	46-54	N1880523
	4	38-42	50-58	N1880524
	5	42-46	54-62	N1880525
AK 26-32cm	6	46-50	58-68	N1880526
	2	30-34	42-50	N1880532
	3	34-38	46-54	N1880533
	4	38-42	50-58	N1880534
AK 26-32cm	5	42-46	54-62	N1880535
	6	46-50	58-68	N1880536

Technical Information

Indications

- Above-knee, below-knee amputations, knee disarticulation
- Sensitive residual limb conditions with a special need for protection
- Phantom pain/sensations/idiopathic residual limb pain
- Reduction of oedema/swelling of the lower limb after amputation

Contraindications

- If the medical findings forbid compression (e.g. heart failure)

Material	Umbrellan® with Clima-Comfort®
Size / Length	See Ordering Information
Amputation Level	Trans tibial and transfemoral
Cleaning	At 30° hand wash
Compression Class	CCL2