

### Skin-Fix Donning Sleeve

Skin-Fix is to be donned on the residual limb before slipping into the socket. Due to the low friction design, donning of the socket becomes an easy and comfortable operation.

Finally, the Skin-Fix can be pulled out through the valve easily.

### Features and Benefits

- Low friction
- Long-life durability
- Simple application

### Ordering Information

Part Number	Size	Upper Circumference	Lower Circumference	Length
NF011749	X-Large	820mm	320mm	1000mm
NF011750	Large	720mm	310mm	1000mm
NF011751	Medium	600mm	260mm	1000mm
NF011752	Small	460mm	260mm	1000mm
NF011753	X-Small	420mm	240mm	1000mm

