

Streifyderm Gel Socks

Streifyderm
 The Streifyderm Gel Sock consists of a nylon fabric sock lined with a 2mm polymer gel, particularly suitable for the moderately active transtibial amputee.

Streifyderm Cushion
 Aimed at meeting the needs of the more active lower extremity amputees, or those with more bony or sensitive residual limbs, the Streifyderm Cushion uses an elasticated nylon sock and a thicker 3mm polymer gel lining.

Streifyderm Endurance
 In order to meet the needs of the very active lower extremity amputee, the Streifyderm Endurance consists of two elasticated nylon socks with the 2mm polymer gel sandwiched between them.

Features and Benefits

- Polymer gel contains a premium medical mineral oil
- Improved skin care as the oil is released
- Reduced skin friction and abrasion
- Increased cushioning and comfort
- Washable
- Suitable for use in many types of socket, regardless of construction or the material used



Streifyderm Gel Socks Ordering Information

Part Number	Size	Length of Gel (cm)	Length of Residual Limb (cm)	Circ Distal (cm)	Circ Proximal (cm)
NMA26A10/1	1	20	25-34	17-24	17-26
NMA26A10/2	2	25	38-44	18-28	18-32
NMA26A10/3	3	25	38-44	20-34	20-38
NMA26A10/4	4	33	43-49	20-34	20-38
NMA26A10/5	5	25	43-49	23-37	23-43
NMA26A10/6	6	33	46-54	23-37	23-43
NMA26A10/7	7	33	46-54	28-42	28-49
NMA26A10/8	8	33	46-54	31-50	31-56

Streifyderm Cushion Ordering Information

Part Number	Size	Length of Gel (cm)	Length of Residual Limb (cm)	Circ Distal (cm)	Circ Proximal (cm)
NMA26A20/1	1	20	25-34	17-24	17-26
NMA26A20/2	2	25	38-44	18-28	18-32
NMA26A20/3	3	25	38-44	20-34	20-38
NMA26A20/4	4	33	43-49	20-34	20-38
NMA26A20/5	5	25	43-49	23-37	23-43
NMA26A20/6	6	33	46-54	23-37	23-43
NMA26A20/7	7	33	46-54	28-42	28-49
NMA26A20/8	8	33	46-54	31-50	31-56

Streifyderm Endurance Ordering Information

Part Number	Size	Length of Gel (cm)	Length of Stump (cm)	Circ Distal (cm)	Circ Proximal (cm)
NMA26A30/1	1	20	25-34	17-24	17-26
NMA26A30/2	2	25	38-44	18-28	18-32
NMA26A30/3	3	25	38-44	20-34	20-38
NMA26A30/4	4	33	43-49	20-34	20-38
NMA26A30/5	5	25	43-49	23-37	23-43
NMA26A30/6	6	33	46-54	23-37	23-43
NMA26A30/7	7	33	46-54	28-42	28-49
NMA26A30/8	8	33	46-54	31-50	31-56

Measure the length of the residual limb from just above the transfemoral condyles, the distal circumference at approximately 4cm from the distal end and the proximal circumference at the patella. Choose the sock appropriate to these measurements, checking that the length of the gel indicated is adequate to allow it to come above the top edge of the socket.