

Wrist Units - Child's 2-Knob Wrist



Child's 2-Knob Wrist Unit

Where regular interchange of a terminal device, rotation, or a lock is required, the Child's 2-Knob Wrist provides a durable unit that can be easily maintained. Commonly used alongside cosmetic and mechanical hands. A comprehensive range of hand and adapter plates are available to allow rotation with a locking feature (see below). Supplied with parallel housing and three screws.

Part Number	Wrist Size	Description	Weight	Height
B05143	40mm (1 1/2")	(Copper) Quick release, lock, 12 position rotation	56g (2 oz)	33mm (1 1/3")
B05142	45mm (1 3/4")	(Black) Quick release, lock, 12 position rotation	81g (3 oz)	33mm (1 1/3")
NSC0375	40mm (1 1/2") + 45mm (1 3/4")	Replacement Screws		



Child's 2-Knob Wrist Hand Plate

The Child's 2-Knob Wrist Hand Plate fits a wide variety of terminal devices to ensure that this wrist is compatible with most fittings. Can be used with cosmetic and mechanical hands.

Part Number	Wrist Size	Description	Weight	Height
B07030	40mm (1 1/2" x 1/2" x 3/8")	1 1/2" Child 2-Knob Rotary	25g (1 oz)	24mm (1")
B07029	45mm (1 3/4" x 1/2" x 3/8")	1 3/4" Child 2-Knob Rotary	31g (1 oz)	24mm (1")



Child's 2-Knob Wrist Adapter Plate

The Child's 2-Knob Wrist Adapter Plate fits a wide variety of terminal devices with a threaded Stem to ensure that this wrist is compatible with most fittings.

Part Number	Wrist Size	Internal Thread	Weight	Height
B00897	40mm (1 1/2" x 1/2" x 3/8")	5/16" BSF (UK Juvenile Standard)	20g (3/4oz)	27mm (1")
B00821	45mm (1 3/4" x 1/2" x 3/8")	5/16" BSF (UK Juvenile Standard)	27g (1 oz)	27mm (1")
B00905	45mm (1 3/4" x 1/2" x 3/8")	3/8" BSW (UK Adult Standard)	43g (1 1/2 oz)	38mm(13/16")