



Seal-In Kit plus Electrode		
Part Number	Description	Power Supply
1	ELSK50	Seal-In Kit Myoelectric 50Hz - Includes 50Hz Electrode
	ELSK60	Seal-In Kit Myoelectric 60Hz - Includes 60Hz Electrode

Check local power supply frequency

Seal-In Kit		
Part Number	Description	Power Supply
2	ELSK	Seal-In Kit

Control Options
 These myoelectrodes offer state-of-the-art electronics for outstanding capture and amplification from signals as low as 10uV. Biocompatible titanium skin contacts provide superior conductivity. Suitable for either Child or Adult arm systems.

Part Number	Description	Power Supply
3	ELEC50	Myoelectric 50Hz
3a	ELEC60	Myoelectric 60Hz

Electrode Cables
 To allow connection between ELEC50 / ELEC60 Electrode and Quick Disconnect Wrist connection or Connector Block Cable. 3-Pin termination to cable.

Part Number	Description	Length
4	ECA-300	Electrode Cable
	ECA-600	Electrode Cable
	ECA-1000	Electrode Cable
4a	B24570	4-Pin Round Elbow Cable - Pink
4b	CBBHA72	Connector Block Cable for Lithium Polymer Battery Receptacle BHA72

Electrode Spare Parts		
Part Number	Description	Power Supply
5	ELECSP	Electrode Forming Kit
5a	B31058	Electrode Suspension Legs
5b	B30819	Electrode Cable Retaining Plug

Silicone Grease		
Part Number	Description	Power Supply
6	NLU0010	Silicone Grease