



Enhanced Steeplon Hand

Strong and robust in construction with four rigid fingers and a spring loaded thumb for manual operation providing a lateral grip.

Part Number	Size	UK Size	Left / Right	Application	Wrist
1 C22047#	7 ¼	2 ¾" x 5 ½" x 2 ⅛"	Left	Adult	54mm (2 ⅛")
C22048#	7 ¼	2 ¾" x 5 ½" x 2 ⅛"	Right	Adult	54mm (2 ⅛")
C22049#	7 ½	3" x 5 ½" x 2 ⅛"	Left	Adult	54mm (2 ⅛")
C22050#	7 ½	3" x 5 ½" x 2 ⅛"	Right	Adult	54mm (2 ⅛")
C22051#	7 ¾	3 ¼" x 5 ½" x 2 ⅛"	Left	Adult	54mm (2 ⅛")
C22052#	7 ¾	3 ¼" x 5 ½" x 2 ⅛"	Right	Adult	54mm (2 ⅛")

Note: Supplied with handplate B01172. Other handplates available on request.

Note: #Denotes finish colour. Add PK for Flesh Pink, ST for Suntan, MB for Mid Brown or DB for Dark Brown.

Colours for Steeplon Hand (for indication only).



Fitting Screw

Part Number	Description	Supplied
2 NSC6011	Replacement Hand Plate Fitting Screw - Steeplon Hand	Each

Leather Gloves

Part Number	Hand Size	UK Glove Size	Colour	Supplied
3 NGL0001	7 ¼ (2 ¾")	6 ½	Brown	Pair
NGL0003	7 ½ (3")	7	Brown	Pair
NGL0005	7 ¾ (3 ¼")	7 ½	Brown	Pair
NGL0006	8 (3 ½")	8	Brown	Pair
NGL0007	8 ¼ (3 ¾")	8 ½	Brown	Pair
NGL0008	9 (4 ½")	9	Brown	Pair
NGL0009	9 ¼ (4 ¾")	9 ½	Brown	Pair
NGL0002	10 (5 ½")	10	Brown	Pair

Note: A bespoke glove service is available for partial hand patients.