


## Wrist Units - Wooden and Fixed Wrists



### Wooden Fixed Wrist

A lightweight wooden wrist that can be modified to produce a bespoke wrist profile. Commonly used alongside a cosmetic hand to create a lightweight prosthesis.

Part Number	Size	Description	Weight	Height
N621020100	60mm (2 3/8")	M12 thread. Parallel, bonded into forearm	40g (1 1/2")	36mm(127/64")



### Fixed Wrist Connector and Former - Round

Commonly used in the manufacture of lightweight cosmetic devices, the Fixed Wrist allows the bonding of a Steeper cosmetic hand to a custom forearm. The wrist former allows the fixed wrist connector to be retrofitted once an outer laminate forearm has been manufactured.

The round fixed wrist connector is also compatible with the Multi-Axial Wrist Housing and Free Flexion Ball Insert (p.109), providing an exoskeletal wrist with variable friction control over rotation, flexion and extension.

Fixed Wrist Part Number	Former Part Number	Wrist Size	Weight	Height
<b>1</b> B05963	<b>2</b> SB12136/1	40mm (1 1/2")	7g (1/4 oz)	16mm (5/8")
B05441	SB12136/2	45mm (1 3/4")	9g (5/16 oz)	19mm (3/4")
B05440	SB12136/3	50mm (2")	11g (3/8 oz)	19mm (3/4")
B05439	SB12136/4	54mm (2 1/8")	13g (1/2 oz)	19mm (3/4")

## Wrist Units - Fixed Wrist and Hand Plates



### Fixed Wrist Connector and Former - Oval

Commonly used in the manufacture of lightweight cosmetic devices, the Fixed Wrist allows the bonding of a Steeper Cosmetic Hand to a custom forearm. The wrist former allows the fixed wrist connector to be retrofitted once an outer laminate forearm has been manufactured.

Fixed Wrist Part Number	Wrist Size	Weight	Height	Former Part Number
<b>3</b> B06252	45mm (1 3/4")	7g (1/4 oz)	22mm (3/4")	<b>4</b> B06986
B06256	50mm (2")	10g (3/8 oz)	22mm (3/4")	B06985
B06255	54mm (2 1/8")	13g (1/2 oz)	32mm (1 1/4")	B06984



### Hand Plate - North American - 1/2"

For use with Foam Hand with Oval base (2 1/2" - 3 3/4").

1. Is used with Steeper foam hand with round interface.
2. Is used with foam hand with oval wrist exposed.

Part Number	Wrist Size	External Thread	Shape	Weight	Height
<b>1</b> B10530	40mm	1/2" x 20 TPI	Round - Steeper Hand Attachment	10g (1/3 oz)	26mm (1")
B22758	40mm	1/2" x 20 TPI	Round - Ottobock Hand Attachment	14g (1/2 oz)	26mm (1")
B22757	45mm	1/2" x 20 TPI	Round - Steeper Hand Attachment	13g (1/2 oz)	26mm (1")
B22756	50mm	1/2" x 20 TPI	Round - Steeper Hand Attachment	16g (1/2 oz)	26mm (1")
B22755	54mm	1/2" x 20 TPI	Round - Steeper Hand Attachment	17g (2/3 oz)	26mm (1")

Both B22757, B22756 and B22755 are also compatible with Ottobock fastening arrangement.

<b>2</b> B06257	45mm	1/2" x 20 TPI	Oval	9g (3/8 oz)	20mm (3/4")
-----------------	------	---------------	------	-------------	-------------



### Hand Plate - European - M12

For use with cosmetic and mechanical hands with a round interface.

Part Number	Wrist Size	External Thread	Shape	Weight	Height
B22768	40mm (1 1/2" x 3/4" x M12)	M12 x 1.5	Round	11g (1/2 oz)	18mm(3/4")
B07786	45mm (1 3/4" x 3/4" x M12)	M12 x 1.5	Round	32g (1 oz)	20mm(3/4")
B06942	50mm (2" x 3/4" x M12)	M12 x 1.5	Round	43g (1 1/2 oz)	20mm(3/4")

Both B06942 and B07786 are also compatible with Ottobock fastening arrangement.