


Infant Foam Filled Cosmetic Gloves



An infant's first cosmetic prosthesis can be fitted between two and six months of age.

This lightweight solution of foam filled, PVC cosmetic glove can assist child development by aiding sitting balance, crawling and play activities.



Infant Foam Filled PVC Gloves

Part Number	Width across Knuckles	Palm Circumference	Middle Finger Length	Application	Left / Right	Length of Foam Fill (Thumb tip to base of glove)
1 AH45/60#	1 1/2"	96mm (3 3/4")	45mm (1 3/4")	Infants	Left	6"
AH45/61#	1 1/2"	96mm (3 3/4")	45mm (1 3/4")	Infants	Right	6"
2 AH45/25#	1 5/8"	98mm (3 7/8")	46mm (1 7/8")	Infants	Left	7"
AH45/26#	1 5/8"	98mm (3 7/8")	46mm (1 7/8")	Infants	Right	7"
3 AH45/64#	1 5/8"	106mm (4 1/8")	35mm (1 3/8")	Infants, Chubby	Left	7"
AH45/65#	1 5/8"	106mm (4 1/8")	35mm (1 3/8")	Infants, Chubby	Right	7"
4 AH45/13#	1 3/4"	105mm (4 1/8")	49mm (2")	Infants	Left	8"
AH45/14#	1 3/4"	105mm (4 1/8")	49mm (2")	Infants	Right	8"
5 AH45/62#	1 3/4"	112mm (4 1/2")	43mm (1 3/4")	Infants, Chubby	Left	8"
AH45/63#	1 3/4"	112mm (4 1/2")	43mm (1 3/4")	Infants, Chubby	Right	8"
6 AH45/27#	1 7/8"	116mm (4 5/8")	54mm (2 1/8")	Infants	Left	8"
AH45/28#	1 7/8"	116mm (4 5/8")	54mm (2 1/8")	Infants	Right	8"
7 AH45/19#	2"	125mm (4 7/8")	57mm (2 1/4")	Infants	Left	8 1/4"
AH45/20#	2"	125mm (4 7/8")	57mm (2 1/4")	Infants	Right	8 1/4"

Note: # Denotes colour. Available in Standard PVC only. 19 colours available - Colour swatch E24654 applies. The above sizes are only available as Foam Filled Gloves.

Manufacturing Guide

