



Friction Housing

Part Number	Description	Thread	Weight	Height
13 B22505	Multi-Axial Wrist Housing 54mm (2 1/8") - Attaches to 54mm Hand Connection Plate. Ball Insert supplied separately		48g (1 1/2 oz)	25mm (1")
14 B22763	Free Flexion Ball Insert - Allows rotation, flexion and extension	Internal Thread 1/2" UNF x 20 TPI	19g (1/2 oz)	24mm (1")
B22503	Free Flexion Ball Insert - Allows rotation, flexion and extension	Internal Thread 3/8" BSW	20g (1/2 oz)	24mm (1")
15 B22764	Locked Flexion Ball Insert - Allows rotation only	Internal Thread 1/2" UNF x 20 TPI	21g (1/2 oz)	26.5mm (1")
B22765	Locked Flexion Ball Insert - Allows rotation only	Internal Thread 3/8" BSW	23g (1/2 oz)	26.5mm (1")

Hand Plates - (North American) For use with Modular Arm System Adapters

Part Number	Wrist Size	External Thread	Shape	Weight	Height
16 B22758	40mm (1 1/2")	1/2" x 20 TPI	Round - Ottobock Hand Attachment	16g (1/2 oz)	25mm (1")
B22757	45mm (1 3/4")	1/2" x 20 TPI	Round - Steeper Hand Attachment	16g (1/2 oz)	25mm (1")
B22756	50mm (2")	1/2" x 20 TPI	Round - Steeper Hand Attachment	17g (1/2 oz)	25mm (1")
B22755	54mm (2 1/8")	1/2" x 20 TPI	Round - Steeper Hand Attachment	19g (1/2 oz)	25mm (1")

B22758, B22757, B22756 and B22755 are also compatible with Ottobock fastening arrangements

17 B06257	45mm (1 3/4") Oval	1/2" x 20 TPI	Oval Wrist Passive Hand	10g (1/2 oz)	25mm (1")
-----------	--------------------	---------------	-------------------------	--------------	-----------

Threaded Insert

Part Number	Description
18 B26431	Threaded Stem. M12 x 1.5mm - 1/2" x 20 TPI
B26430	Threaded Stem. M12 x 1.5mm - M12 x 1.5mm

Inner Foam Cosmesis

Part Number	Description	Weight
19 C11151	Inner Foam Cosmesis - Transhumeral	134g (5 oz)
C11152	Inner Foam Cosmesis - Transradial	134g (5 oz)

Outer Fabric Cosmesis

Part Number	Description	Thread	Length
20 NOB0002	Stockinette Cover	Flesh Pink	1m lengths