



Conversion Sleeve		
Part Number	Description	Supplied
1 B06005	Conversion Sleeve - External Thread 1/2" UNF x 20 TPI, Internal Thread 3/8" BSW To convert UK External Thread to US External Thread	Each
B21871	Conversion Sleeve - External Thread M12 x 1.5mm, Internal Thread 3/8" BSW To convert UK External Thread to European External Thread	Each
B23372	Conversion Sleeve - External Thread 1/2" UNF x 20 TPI, Internal Thread 5/8" BSF To convert Child UK External Thread to US External Thread	Each


Light Tool Holder					
An adjustable clamp with two jaws for gripping around handled tools. The head may be set at various angles and jaws tightened by means of a wing nut.					
Part Number	Description	Left / Right	Thread	Length	Weight
2 B05191	Light Tool Holder - Adult	Left or Right	3/8" BSW	120mm (4 3/4")	250g (9 oz)

Spade Grip					
Suitable for use with gardening or sports equipment. Provides a full range of unrestricted wrist movement.					
Part Number	Description	Left / Right	Thread	Length	Weight
3 C07077	Spade Grip - Adult	Left or Right	3/8" BSW	125mm (5")	362g (12 1/2 oz)

Universal Joint					
A quick release 2 1/2" adapter which provides unrestricted wrist rotation with no lock and free movement in all directions.					
Part Number	Description	Left / Right	Thread	Length	Weight
4 B00774	Universal Joint - Adult	Left or Right	3/8" BSW	70mm (2 3/4")	131g (4 1/2 oz)
4a SB18012	Universal Joint - Adult	Left or Right	3/8" BSW Internal and External		131g (4 1/2 oz)

Supplied with 3/8" BSW External Thread as opposed to adapter. This allows choice in selection of interface. Provides unrestricted wrist rotation with no lock, and free movement in all directions.

**Note:** A full range of adapter plates are available. Please see page 111. Thread locking compound must be used when assembling adapters to functional devices.



**Fishing Rod Holder**

Used to support a fishing rod, allowing the sound hand freedom to control the reel.

Part Number	Description	Left / Right	Thread	Length	Weight	
5	B05882	Fishing Rod Holder	Left or Right	3/8" BSW	140mm (5 1/2")	143g (5 oz)

**Walking Stick Holder**

The walking stick holder is connected to the prosthesis with an angled tube. Supplied with 2 1/2" adapter plate. Walking Stick not included.

Part Number	Description	Left / Right	Thread	Length	Weight	
6	B05640	Walking Stick Holder - Adult	Left or Right	3/8" BSW	140mm (5 1/2")	200g (7 oz)

**Typing Appliance**

Supplied with an angled head and rubber tip for operation of a keyboard.

Part Number	Description	Left / Right	Thread	Length	Weight	
7	B05362	Typing Appliance - Adult	Left or Right	3/8" BSW	115mm (4 1/2")	48g (1 1/2 oz)
	B05481	Replacement Rubber Tip	n/a	n/a	n/a	n/a

**Plain Hook**

For use in manual handling and lifting applications.

Part Number	Description	Left / Right	Thread	Length	Weight	
8	B05138	Plain Hook - Adult Short	Left or Right	3/8" BSW	135mm (5 3/8")	109g (4 oz)
	B05140	Plain Hook - Adult Long	Left or Right	3/8" BSW	300mm (11 3/4")	169g (6 oz)

**'C' Hook**

A multi-purpose hook for use with gardening and farming tools (e.g. broom, shovel, fork, wheelbarrow) where pushing and pulling are required. The head can swivel and flex, providing a more adaptive device. The plain 'C' Hook retains functional head shape but is attached to the stem in a fixed position.

Part Number	Description	Left / Right	Thread	Length	Weight	
9	B00551	'C' Hook - Adult	Left or Right	3/8" BSW	125mm (5")	260g (9 oz)
	B00770	Plain 'C' Hook - Adult	Left or Right	3/8" BSW	140mm (5 1/2")	240g (8 1/2 oz)

**Nail Brush Holder**

This stainless steel adapter has two drilled holes for attachment to the object of a user's choice, e.g. nailbrush, wire brush, files, paint roller etc. Supplied as shown.

Part Number	Description	Left / Right	Thread	Length	Weight	
10	B05353	Nail Brush Holder - Adult	Left or Right	3/8" BSW	100mm (4")	52g (2 oz)

**Note:** A full range of adapter plates are available. Please see page 111. Thread locking compound must be used when assembling adapters to functional devices.



**Ball Pein Hammer**  
Suitable for light metal work and general engineering use.

Part Number	Description	Left / Right	Thread	Length	Weight
11	B05210	Ball Pein Hammer - Adult	3/8" BSW	130mm (5 1/8")	342g (12 oz)

**Claw Hammer**  
Suitable for woodwork and general use.

Part Number	Description	Left / Right	Thread	Length	Weight
12	B05212	Claw Hammer - Adult	3/8" BSW	130mm (5 1/8")	599g (1lb 5 oz)

**Potato Holder**  
Allows items such as vegetables, fish, bread and meat etc. to be held securely during food preparation, while leaving the sound hand free to peel or cut.

Part Number	Description	Left / Right	Thread	Length	Weight
13	B00538	Potato Holder - Adult	3/8" BSW	120mm (4 3/4")	97g (3 1/2 oz)

**Table Spoon**  
Large standard spoon manufactured in stainless steel.

Part Number	Description	Left / Right	Thread	Length	Weight
14	B05356	Table Spoon - Adult	3/8" BSW	140mm (5 1/2")	42g (1 1/2 oz)
14a	NGO3F82L	Angled Table Spoon - Adult	3/8" BSW	175mm (6 7/8")	82g (3 oz)
	NGO3F82R	Angled Table Spoon - Adult	3/8" BSW	175mm (6 7/8")	82g (3 oz)

**Table Knife**  
Large standard knife manufactured in stainless steel.

Part Number	Description	Left / Right	Thread	Length	Weight
15	B05354	Table Knife - Adult	3/8" BSW	160mm (6 1/4")	37g (1 1/2 oz)
15a	NGO3F81	Steak Knife - Adult	3/8" BSW	130mm (5 1/8")	40g (1 1/2 oz)

**Note:** A full range of adapter plates are available. Please see page 111.  
Thread locking compound must be used when assembling adapters to functional devices.



**Table Fork**  
Large standard fork manufactured in stainless steel.

Part Number	Description	Left / Right	Thread	Length	Weight
16 B05355	Table Fork - Adult	Left or Right	3/8" BSW	140mm (5 1/2")	40g (1 1/2 oz)
16a NGO3F83L	Angled Table Fork - Adult	Left	3/8" BSW	170mm (6 3/4")	64g (2 1/4 oz)
NGO3F83R	Angled Table Fork - Adult	Right	3/8" BSW	170mm (6 3/4")	64g (2 1/4 oz)

**Vegetable Peeler**  
Sharp cutting edge allows easy peeling of fruit or vegetables.

Part Number	Description	Left / Right	Thread	Length	Weight
17 B11476	Vegetable Peeler - Adult	Left or Right	3/8" BSW	110mm (4 3/8")	16g (1/2 oz)

**Cue Holder**  
Lined in rubber with an adjustable clamp to securely grip the cue.

Part Number	Description	Left / Right	Thread	Length	Weight
18 B05700	Cue Holder - Adult	Left or Right	3/8" BSW	90mm (3 1/2")	16g (5 1/2 oz)


**Cue Rest**  
Designed to support the cue, this device has a rotating head to provide a stable base and three different cueing heights.

Part Number	Description	Left / Right	Thread	Length	Weight
19 B07318	Cue Rest - Adult	Left or Right	3/8" BSW	165mm (6 1/2")	100g (3 1/2 oz)

**Telescopic Dressing Stick**  
Extendable device designed to facilitate buttoning and unzipping using either mouth or residual digits. May be used with mouthpiece.

Part Number	Description	Left / Right	Weight
20 B05553	Telescopic Dressing Appliance	n/a	75g (2 1/2 oz)
20a SB12202	Mouthpiece	n/a	15g (1/2 oz)

**Note:** A full range of adapter plates are available. Please see page 111.  
Thread locking compound must be used when assembling adapters to functional devices.



**Motoring Ball and Clamp - Adult**

For driving applications. A rim-mounted single clamp arrangement fits most types of steering wheels without compromising deployment of safety devices.

Part Number	Description	Left / Right	Note	Weight
21 E22252	Motoring Ball and Clamp - Adult	n/a	Available UK only	382g (13 ½ oz)
E22252A	Motoring Clamp with Link - Adult	n/a	Available UK only	255g (9 oz)

**Note:** E22252A is supplied with link bar to allow connection to a second Motoring Ball and Clamp mounted on the steering wheel.

**Ball and Clamp**

This device may be used as a handlebar-mounted aid for adults and children wishing to cycle.

Part Number	Description	Left / Right	Note	Weight
21a SP118	Ball and Clamp	n/a	Not for use with motorised vehicles	115g (4 oz)

**Driving Cup**

For driving applications. This device is screwed into a palm mounted fitting. It must be used with an appropriate steering ball and clamp. A fitted PVC cover gives a secure level of attachment while allowing rapid removal if necessary.

Part Number	Description	Left / Right	Note	Thread	Weight
22 B00447	Driving Cup - Adult	n/a	Available UK only	¾" BSW	126g (4 ½ oz)

**Note:** Alternative cup sizes may be manufactured. Please contact customer services with your requirements.

**Driving Cup on Stem**

For driving applications. Must be used with an appropriate steering ball and clamp. A fitted PVC cover gives a secure level of attachment while allowing rapid removal if necessary.

Part Number	Description	Left / Right	Note	Thread	Weight
23 B00451	Driving Cup on Stem - Adult	n/a	Available UK only	¾" BSW	126g (4 ½ oz)

**Note:** Alternative cup sizes may be manufactured. Please contact customer services with your requirements.

**Extension Bar**

An extension piece which may be used to increase length between two prosthetic components.

Part Number	Description	Left / Right	Thread	Weight
24 B00445	Extension Bar 3" - Adult	n/a	¾" BSW External and Internal	71.5g (2 ½ oz)
B00446	Extension Bar 4" - Adult	n/a	¾" BSW External and Internal	95.5g (3 ½ oz)

**Note:** A full range of adapter plates are available. Please see page 111. Thread locking compound must be used when assembling adapters to functional devices.