



# TASKA CX™

Faster. Stronger. More connected.

 steeper  **TASKA™**

by eqwal.

# TASKA's most advanced prosthetic hand



## 23 Grips

Complete most tasks with 3 core grips



## 20kg

Overall carry load



## 0 - 110mm

Maximum opening range  
\*dependent upon hand size



## IP67

Waterproof rating

### Ground Breaking Motor Design

Perform a wider range of tasks faster, with more precision and superior grip strength. Whether it's shaking hands, holding a broom or using utensils

### Capacitive Touch

All fingertips, including the thumb, have capacitive touch and can be used on any common touch screen

### One-click Grip Tip Replacement

Grip Tips can be easily replaced with a single simple click to keep you moving quickly

### Intuitive Button Layout

Easily cycle through grips, control Bluetooth® and more with the ergonomically designed button panel

### Coverplate Customisation

Customise the look of your device with replaceable Coverplates available in a range of colours and designs

### Anatomical Design

Contoured lines and discreet Wrist Flexion buttons create a streamlined profile, ideal for those seeking a sleek design




*This new hand has  
been a real jump in  
quality and trust for  
me.*

Korrin Barrett



# Complete almost any task

With 18 pre-programmed grips and the ability to create five of your own, completing almost any task is easy. The grips you will use for most of your daily activities are based on three fundamental grasps, whilst dynamic, static, and speciality grips are also available.



## Customise your grips using the MyTASKA App

With the MyTASKA app, you're in complete control, with easy access to customisable settings and grips for your hand. Use your hand your way by adding, removing and reordering your grips, creating a configuration that works best for you.



### Customisable Grips


The app features a library of 18 preconfigured grips with the opportunity to add a further five custom options. From precision grips to power grips, you can simply select your grip of choice and customise it to your hand.



### Built-in Games

You can hone signal accuracy and improve control of your TASKA hand with three fun games. Refine clear open and close signals and practice EMG (electromyography) proportional control.

Download the MyTASKA app from your App Store





## No glove, no problem

With an IP67 waterproof rating, you can wash dishes, water the garden or simply walk in the rain, safe in the knowledge you're not putting your hand at risk.

Cleaning your prosthesis with both hands together is easy and means hygiene is significantly improved.

All TASKA hands come with wrist options that include TASKA's HydroSeal™ Waterproofing Technology, allowing submersion past the wrist.

Please note, the hand is only waterproof to the wrist if used with a seal ring. Any components above the wrist are not protected.

## Make it your own with 5 different coverplates

Customise your CX with a choice of two base colours and five interchangeable coverplates. This allows you to express your personality through your prostheses and style it your way.



black

sand



black

grey

dusk

white

sand

Compatible with a wide variety of myoelectric components including:



# Espire™ Elbow

Pro & Hybrid



For a full list of compatible components, visit [www.taskaprosthetics.com/support/compatibility](http://www.taskaprosthetics.com/support/compatibility)

Espire Elbow pictured above with TASKA HandGen2.

Steeper Group  
Unit 3 Stourton Link  
Intermezzo Drive  
Leeds, LS10 1DF  
Email: [customerservices@steepergroup.com](mailto:customerservices@steepergroup.com)  
Tel: +44 (0) 870 240 4133

[www.steepergroup.com](http://www.steepergroup.com)

STPPR334

 **steeper**  
by eqwal.