



Titanium (5XTi)					
Part Number	Description	Left / Right	Thread	Length	Weight
1 NHD0003T	Large Titanium Hook	Left	3/8" BSW	124mm (4 7/8")	151g (5 oz)
NHD0004T	Large Titanium Hook	Right	3/8" BSW	124mm (4 7/8")	151g (5 oz)

Stainless Steel (5X)					
Part Number	Description	Left / Right	Thread	Length	Weight
1 NHD0003	Large Stainless Steel Hook	Left	3/8" BSW	124mm (4 7/8")	213g (7 1/2 oz)
NHD0004	Large Stainless Steel Hook	Right	3/8" BSW	124mm (4 7/8")	213g (7 1/2 oz)

Aluminium (5XA)					
Part Number	Description	Left / Right	Thread	Length	Weight
1 NHD0001	Large Aluminium Hook	Left	3/8" BSW	124mm (4 7/8")	113g (4 oz)
NHD0002	Large Aluminium Hook	Right	3/8" BSW	124mm (4 7/8")	113g (4 oz)

Aluminium (99X)					
Part Number	Description	Left / Right	Thread	Length	Weight
2 NHD0005	Medium Aluminium Hook	Left	3/8" BSW	98mm (3 7/8")	85g (3 oz)
NHD0006	Medium Aluminium Hook	Right	3/8" BSW	98mm (3 7/8")	85g (3 oz)

Aluminium (10X)					
Part Number	Description	Left / Right	Thread	Length	Weight
3 NHD0007	Small Aluminium Hook	Left	5/16" BSF	86mm (3 3/8")	85g (3 oz)
NHD0008	Small Aluminium Hook	Right	5/16" BSF	86mm (3 3/8")	85g (3 oz)

Aluminium (Small)					
Part Number	Description	Left / Right	Thread	Length	Weight
4 B05157	Small Aluminium Hook	Left	5/16" BSF	76mm (3")	42g (1 1/2 oz)
B05158	Small Aluminium Hook	Right	5/16" BSF	76mm (3")	42g (1 1/2 oz)

Note: Adult hooks are supplied with standard UK operating cord connection ring.

Note: For hooks with bonded jaw linings a repair service is available. Order by adding EX to part number e.g. NHD0003EX. Contact your local supplier for return details.

Note: A full range of adapter plates are available. Please see page 111.

Spare Parts		
Part Number	Description	Supplied
5 B11470	Rubber Rings for Hook Tension	Pack of Ten
5a B09439	Large Rubber Rings	Each
6 NRB0050	1/2" Rubber Tube Lining	Per Metre

Jaw Covers for Child Hook, B05157 and B05158		
Part Number	Description	Supplied
7 B19464	Half Round Jaw Cover	Each
E19720	Half Round Jaw Cover	Pack of Twenty
7a B19572	Full Round Jaw Cover	Each
E19927	Full Round Jaw Cover	Pack of Twenty