

Infant & Paediatric Hands

Greek Series

Durable infant hand constructed using soft polymer in a unisex style. They have realistic finger and hand features, and come in three models to reflect the development of a child's hand. The design and shape of the hand means they can be used to hold very light objects, such as plastic toys or dummies, and assist young infants as they become mobile.



| Part Number | Model |
|-------------------|-------|
| NALPHA- [] - [] | Alpha |
| NBETA- [] - [] | Beta |
| NGAMMA- [] - [] | Gamma |

[] = LEFT or RIGHT

[] = Colour (Tan = T, Medium Brown = MB, Dark Complexion = DK)

Features

- Easy to clean
- Elastic, flexible holding capability with no cable required
- ½" x 20TPI removable threaded stud

Application

- All levels of upper limb absence
- Ages 4 months to 3 years old

Technical Information

| Model | Alpha | Beta | Gamma |
|---------------|-------------|--------------|--------------|
| Application | 4-12 months | 10-24 months | 20-36 months |
| Length | 5cm | 6.4cm | 7cm |
| Width | 4cm | 4.5cm | 4.8cm |
| Grasp-Opening | 3.8cm | 4.3cm | 4.6cm |
| Base Diameter | 3.3cm | 3.7cm | 3.7cm |
| Weight | 40g | 45g | 48g |

Infant & Paediatric Hands

Infant 2 Hand



The Infant 2 is a passive prosthetic hand, suitable for infants beginning to crawl or walk. The unique 'cup' shape of the palmar surface aids in pulling and pushing objects, whilst the flexible, durable polymer material is ideal for bi-lateral applications as it enables a grasping action.



| Part Number |
|---------------------|
| NINFANT2- [] - [] |

[] = Side (L = Left, R = Right)

[] = Colour (Tan = T, Medium Brown = MB, Dark Complexion = DK)

Features

- No cable required
- Soft material for ambulating infants
- Easy to clean
- Realistic finger and thumb design features
- ½" x 20TPI (M12 available, order conversion part number: NCUSTM)

Application

- Infants aged 6-18 months
- All levels of upper limb absence

Technical Information

| | |
|---------------------|-------|
| Length | 6.6cm |
| Width | 5cm |
| Weight | 57g |
| Wrist Diameter | 3.8cm |
| Wrist Circumference | 10cm |