

Downhill Racer Ski



A skiing device for snapping-on to ski poles, which features a pre-flexed angle for snow clearance and optimised control. It is constructed using flexible, high-performance polyurethane for durability, and comes in multiple flexibilities to meet the user's needs.

Part Number

NSKI-RACER- - *

= Side (LEFT or RIGHT)
 = Colour (Black = *leave blank*, TL = Teal, BLU = Blue, V = Violet)

*Please specify age, weight and performance criteria when ordering as these devices are built to individual requirements.

Features

- Comes with a press-on rubber pole cap for protection
- Designed to fit ~1.6cm diameter ski poles, and can also be used on hiking poles
- ½" x 20TPI (M12 available, order conversion part number: NCUSTM)

Applications

- Youth and adult
- Suitable for downhill and cross-country skiing

Technical Information

Length	8.4cm
Width	6.4cm
Thickness	3.3cm
Weight	114g

Ski 2



The Ski 2 is a terminal device which can be used with or without a cable for quick and precise pole plants, and utilises a quick disconnect pole system. It has been designed to assist with balance and body posture when skiing to ensure users are safe.

This device is optimal for nordic skiing where upper body power is required for propelling the body forward.

Part Number

NSKI2TD**

Technical Information

Length	8.4cm
Width	3.8cm
Weight	298g
Materials	Aluminium, stainless steel, rubber, plastic

**Please specify whether user is a novice or advanced skier when ordering as two types of shock absorbing mounts are available.



Features

- ½" x 20TPI (M12 available, order conversion part number: NCUSTM)
- Includes a shock absorbing flexion mount
- Accepts 0.5cm and 0.7cm diameter cable ball fittings
- Can be used with poles with a diameter from 1.5 - 1.9cm

Application

- Cross-country (nordic) and downhill (alpine) skiers, from novice to advanced
- For youth and adults

Rocky Mountain High Skip Ropes



The Rocky Mountain High Skip Ropes are set of lightweight, high-performance skipping ropes along with an optimally angled adapter which screws into a standard or disconnect-style wrist adapter. The rope length can be adjusted easily to suit individual requirements, and fits securely onto the rope adapter.

Part Number	Model
NJUMP-ROPE	Skip Ropes
NJR-ADAPTER	Skip Ropes Adapter

Features

- ½" x 20TPI adapter (M12 available, order conversion part number: NCUSTM)

Application

- Children aged 8 up to adults
- Right or left application - Bilateral custom sets can be made (please contact Customer Services to enquire)

Technical Information

Length	13cm
Weight	128g
Colour	Aluminium, black & white
Materials	Nylon, neoprene, aluminium, steel